




Peak Flow Meters



Patients can be confident in the proven, time tested performance of the Mini-Wright™ and Mini-Wright AFS™ peak flow meters, or the economical, accurate Airzone™ peak flow meter.

Years of clinical tests have proven that both the Standard Range Mini-Wright™ and Low Range Mini-Wright AFS™ peak flow meters have delivered exceptional and consistent accuracy.

The Airzone™ peak flow meter, provides patients an economical choice of accurate peak flow monitoring.



Manage your asthma

- All meters for use by children and adults
- All individually calibrated
- Mini-Wright models have built-in, one-way valves
- Optional Windmill Trainers for Mini-Wright models, teach proper peak flow technique and improve monitoring compliance. Ideal for children as young as 2-1/2 years
- Airzone folding handle adds to portability and convenience
- Airzone's tape zone markers prevent inadvertent repositioning

Choose accurate peak flow monitoring as part of your asthma management program.

Peak Flow Meters

Mini-Wright™ & Mini-Wright AFS™ Specifications

| | |
|----------------------------|--|
| ATS Scale | |
| MINI-WRIGHT Standard | 60-880 L/min |
| MINI-WRIGHT AFS Low | 30-400 L/min |
| Scale Increments | |
| | 10LPM |
| Size | |
| MINI-WRIGHT Standard | 5-7/8"H X 1-7/8"W |
| MINI-WRIGHT AFS Low | 5-7/8"H X 1-1/2"W |
| Weight | |
| MINI-WRIGHT Standard | 2.9 oz |
| MINI-WRIGHT AFS Low | 2.2 oz |
| Materials | |
| | All Latex Free |
| Patient Information | |
| | Operating Instructions/ Patient Diary in English, Spanish and French |



PEAK FLOW METER



Ordering Information

| | Order # | Description/Quantity | NHRIC# |
|-----------------------------------|---------|---|----------------|
| Mini-Wright | 3103104 | Standard - each or 10/case | 3 8462103104 0 |
| | 3104001 | Low Range - each or 10/case | 3 8462104001 1 |
| Windmill Trainer | 3122015 | Blue/Red trainer - 5/bag | 3 8462122015 4 |
| Airzone | 3130045 | each or 10/case | 3 8462130045 0 |
| | 3122200 | Adult - 100/bag | 3 8462122200 4 |
| | 3123200 | Pediatric - 100/bag | 3 8462123200 3 |
| Disposable Mouthpieces | 3122064 | One-way exp. - 200/bx | 3 8462122064 2 |
| | 3122120 | Pediatric/Small Mouth Adult - bag of 20 | 3 8462122120 5 |
| Disposable Adapters | 3122220 | Pediatric/Small Mouth Adult - box of 100 | 3 8462122220 2 |



Alliance Tech Medical, Inc.™

PO Box 6024 / Granbury, TX 76049
1.800.848.8923 / 817.326.6357
Fax: 817.326.2182 / www.alliancetechmedical.com

Eastern Office 302.838.2638 **Central Office** 630.910.0993
Southern Office 817.326.3183 **Western Office** 800.848.8923

Airzone™ Specifications

| | |
|----------------------------|---|
| ATS Scale | 60-720 L/min |
| Scale Increments | 10LPM |
| Size | 5-1/2"H X 1-3/4"W |
| Weight | 1.5 oz |
| Materials | All Latex Free |
| Patient Information | Operating Instructions/ Patient Diary in English, Spanish and French |



Disposable Mouthpieces

Economical and convenient. Reduce risk of cross-infection, helping to protect patient and equipment of contamination. One-way valve mouthpieces available for expiratory use.



Mini-Wright™ Windmill Trainer™

The Windmill Trainer is a revolution in teaching proper peak flow technique. By improving technique, use of the Windmill Trainer reduces variability in peak flow readings and improves monitoring compliance. Ideal for use with children as young as 2-1/2 years old, the elderly and those with language or learning difficulties. The Windmill Trainer is an effective yet enjoyable tool. Use the Windmill Trainer with both the Mini-Wright AFS Low Range and Mini-Wright Standard Range peak flow meters.

Call today for Alliance Tech Medical's complete line of COPD, respiratory, asthma management products, sanitizers, splash shields, ventlab chambers and many other items.

Distributed By: